

StemAdhere™**User Protocol****StemAdhere™ Pluripotency Substrate (Product No. S2071-500UG)**

- 500 µg StemAdhere Pluripotency Substrate, sterile filtered, 250 µg/ml
Stability and storage: 6 months at -80 °C or 2 weeks at 4 °C, undiluted

StemAdhere™ Pluripotent Cell Growth Kit Components (Product No. S2112)

- 500 µg StemAdhere Pluripotency Substrate, sterile filtered, 250 µg/ml
Stability and storage: 6 months at -80 °C or 2 weeks at 4 °C, undiluted
- 100 ml Dilution Buffer, sterile filtered
Stability and storage: 6 months at 15-30 °C. Do not refrigerate.
- 100 ml Cell Release Buffer, sterile filtered
Stability and storage: 6 months at 15-30 °C
- 8 ea 6-well polystyrene plates, non-TC treated, sterile individually wrapped




Fig. 1. Colony of H1 cells growing on StemAdhere™ coated non-TC treated 24-well plate

StemAdhere™ Pluripotency Substrate is a defined single-component cell attachment factor intended to support proliferation of human stem cells (ES and iPS) and maintenance of the pluripotent state [1, 2]. StemAdhere consists entirely of human sequences, and is to be used for coating of sterile, non-tissue culture treated polystyrene dishes and plates. Coating should be performed on the day of use; plates with coating solution may also be stored overnight at 4 °C. Cells grown on StemAdhere exhibit a flatter and larger cell morphology than seen on commonly used matrices. This morphology is indicative of greater cell adhesion to the substrate.

Important!

- Coat only NON-tissue culture treated polystyrene plasticware
- Use polypropylene tubes for dilution of StemAdhere
- Use a 5 ml glass pipette when passaging the cells

Primorigen Biosciences • 510 Charmany Drive • Madison, WI 53719

☎ 608.441.8332 • ✉ info@primorigen.com • www.primorigen.com

Trademark and Patent Information: Primorigen®, and Primorigen Biosciences® are registered trademarks of Primorigen Biosciences.
StemAdhere™ is a trademark of Primorigen Biosciences (Patent Pending)

Protocols for use of StemAdhere™

Materials required for coating:

- Non-tissue culture treated polystyrene plasticware (supplied 6-well plates or see page 5 for suitable alternatives)
- Sterile polypropylene conical tubes, 15 or 50 ml
- Sterile transfer pipettes
- Dilution Buffer (Product No. S2072-100ML)

Materials required for establishment and passaging of cells:

- Prewarmed culture medium (e.g. mTeSR[®]1, TeSR™2, or MEF-conditioned ES medium)
- Cell Release Buffer (Product No. S2117-200ML; acceptable results may also be obtained with Lonza Versene #17-711E)
- 5 ml glass pipette, Fisherbrand #13-678-27E

Coating plasticware with StemAdhere

*Note: This protocol should be performed under aseptic conditions in a Biological Safety Cabinet. **USE ONLY NON-TISSUE CULTURE TREATED PLASTICWARE.***

1. Thaw the vial of StemAdhere at room temperature.
2. Dilute StemAdhere to 10 µg/ml with Dilution Buffer in a polypropylene conical tube. Prepare a sufficient volume for the desired plasticware (see Table 1 on page 4).
The required dilution is 1:25, e.g. 240 µl StemAdhere diluted into 5.76 ml Dilution Buffer is sufficient for one 6-well plate. Optimal concentration of StemAdhere may vary depending upon cell type- try 15 and 20 µg/ml if attachment seems weak. Aliquot unused thawed StemAdhere in sterile tubes and refreeze at -80 °C for long term storage.
3. Mix gently (do not vortex).
4. Deposit required volume of diluted StemAdhere in a pool in the center of the well (e.g. 1 ml/well in a 6-well plate). After adding StemAdhere to the desired wells, gently shake the plate horizontally, side to side and forward-backward to spread the coating solution across the entire well or plate surface. Smaller volumes than recommended may be used, but will require more vigorous shaking to completely cover the surface.
5. Incubate 1 hr at 37 °C in a CO₂ incubator.
After 37 °C incubation, plates may be stored overnight at 4 °C (wrap plates in Parafilm[®] or place in ziploc bag to minimize evaporation).
6. Remove StemAdhere solution and rinse once with Dilution Buffer before addition of cells. Do not allow wells to dry out.
7. Medium may be added and the plate held at 37 °C in a CO₂ incubator until addition of the cells (not more than 1 hr).

Primorigen Biosciences • 510 Charmany Drive • Madison, WI 53719

☎ 608.441.8332 • ✉ info@primorigen.com • www.primorigen.com

Trademark and Patent Information: Primorigen®, and Primorigen Biosciences® are registered trademarks of Primorigen Biosciences. StemAdhere™ is a trademark of Primorigen Biosciences (Patent Pending)

Converting cells from Matrigel to growth on StemAdhere™

It is essential to prevent proteolysis of cell surface proteins when preparing colonies of pluripotent cells for establishment on StemAdhere coated plates. The following non-enzymatic steps should be used for colonies growing on Matrigel or other feeder-free ECM substrates. Suggested volumes are for 6-well plates; adjust as needed.

1. Aspirate medium, rinse with 1 ml/well Cell Release Buffer and aspirate.
2. Add 1 ml/well Cell Release Buffer.
3. Incubate at room temperature for **9 minutes**.
4. Gently aspirate Cell Release Buffer.
5. Quickly remove the colonies from each well with 2 ml/well growth medium using a 5 ml glass pipette. Forcefully dispense the entire volume of medium with the pipette at a 90° angle to the plate, starting at the bottom of the well and moving upwards. The goal is to wash the entire content of the wells off the plate with one dispensing (or at most two) of the medium per well. **Avoid scraping**. The colonies should be easily removed without having to scrape the bottom of the plate. **Do not further break up the colonies**
6. Dropwise evenly distribute the colony suspension into each well of a medium containing prepared StemAdhere 6-well plate at desired split ratio (1:3 to 1:6 split recommended, depending on specific cell line)
Optional: ROCK inhibitor and Thiazovivin have been reported to increase the efficiency of pluripotent cell survival after replating [3, 4]. ROCK inhibitor Y-27632 added at 10 μ M final concentration, and Thiazovivin at 2 μ M final concentration may enhance survival of cells and adaptation to StemAdhere.
7. Ensure that each well has at minimum 2 ml/well medium; add additional medium if needed.
8. Refeed colonies daily and monitor morphology and confluency. Passage when adjacent colonies begin to touch.

Converting feeder-dependent cells to growth on StemAdhere™

It is essential to prevent proteolysis of cell surface proteins when preparing feeder-dependent colonies of pluripotent cells for establishment on StemAdhere coated plates. The following manual passaging steps should be used for colonies growing on feeder cells.

1. Aspirate medium and replace with fresh 37 °C growth medium.
2. Manually cut each colony into 4-6 pieces using a sterile 26 g needle (attach to 1 ml syringe for easier handling). Alternatively, use a roller tool to cut the colonies, such as the StemPro® EZPassage™ disposable tool from Invitrogen.
3. Dislodge the colony fragments with a sterile pipette tip and aspirate the floating colonies with the culture medium.
4. Transfer to StemAdhere-coated well or dish of the same size. Do not dilute cells at this first passage (split ratio 1:1).

Optional: ROCK inhibitor and Thiazovivin have been reported to increase the efficiency of pluripotent cell survival after replating [3, 4]. ROCK inhibitor Y-27632 added at 10 μ M final concentration, and Thiazovivin at 2 μ M final concentration may enhance survival of cells and adaptation to StemAdhere.

Primorigen Biosciences • 510 Charmany Drive • Madison, WI 53719

☎ 608.441.8332 • ✉ info@primorigen.com • www.primorigen.com

Trademark and Patent Information: Primorigen®, and Primorigen Biosciences® are registered trademarks of Primorigen Biosciences. StemAdhere™ is a trademark of Primorigen Biosciences (Patent Pending)

5. Refeed colonies daily and monitor morphology and confluency. Passage when adjacent colonies begin to touch.

Passaging cells grown on StemAdhere™

The split ratios of 1:3 to 1:6 from StemAdhere to StemAdhere have resulted in desirable attachment and confluence rates. It is recommended to use these ratios depending on the confluence of the starting wells. Suggested volumes are for 6-well plate, adjust as needed.

1. Aspirate medium from wells, rinse with 1ml/well Cell Release Buffer and aspirate.
2. Add 1 ml/well Cell Release Buffer.
3. Incubate at room temperature for **2-3 minutes**.
4. Gently aspirate Cell Release Buffer.
9. Quickly remove the colonies from each well with 2 ml/well growth medium using a 5 ml glass pipette. Forcefully dispense the entire volume of medium with the pipette at a 90° angle to the plate, starting at the bottom of the well and moving upwards. The goal is to wash the entire content of the wells off the plate with one dispensing (or at most two) of the medium per well. **Avoid scraping.** The colonies should be easily removed without having to scrape the bottom of the plate. **Do not further break up the colonies**
5. Dropwise evenly distribute the colony suspension into each well of a medium containing prepared StemAdhere 6-well plate at desired split ratio (1:3 to 1:6 split recommended, depending on specific cell line)
Optional: ROCK inhibitor and Thiazovivin have been reported to increase the efficiency of pluripotent cell survival after replating [3, 4]. ROCK inhibitor Y-27632 added at 10 μM final concentration, and Thiazovivin at 2 μM final concentration may enhance survival of cells and adaptation to StemAdhere.
6. Ensure that each well has at minimum 2 ml/well medium; add additional medium if needed.
7. Refeed colonies daily and monitor morphology and confluency. Passage when adjacent colonies begin to touch.

Table 1. Recommended volumes for coating

	Wells/dish				Dish/flask size			
	96	24	12	6	60 mm	100 mm	T25	T75
area per well (cm ²)	0.3	2	3.8	9.6	28.3	78.5	25	75
recommended volume per well (ml)	0.05	0.25	0.50	1.0	3.0	6.0	3.0	5.0
minimum volume per well (ml) at 78 μ l/cm ²	0.02	0.16	0.30	0.75	2.20	6.12	1.95	5.85

Table 2. Polystyrene (non-TC treated) plasticware sources

Manufacturer	Type	Product No.	StemAdhere tested
Falcon/BD	24 well polystyrene	351147	Yes
Fisher	60 mm dish	08-757-BA	Yes
Greiner BioOne	6 well polystyrene	657185	Yes
	60 mm dish	628162	Yes

Available separately

StemAdhere Pluripotency Substrate	S2071-1MG
Cell Release Buffer (100 ml)	S2072-100ML
Dilution Buffer (100 ml)	S2117-100ML
6-well Polystyrene Plates (100/cs)	S2118-1CS

References

1. Nagaoka, M., K. Si-Tayeb, T. Akaike, and S.A. Duncan, *Culture of human pluripotent stem cells using completely defined conditions on a recombinant E-cadherin substratum*. BMC Dev Biol. 10(1): p. 60.
2. Chen, T., D. Yuan, B. Wei, J. Jiang, J. Kang, K. Ling, Y. Gu, J. Li, L. Xiao, and G. Pei, *E-cadherin-Mediated Cell-Cell Contact is Critical for Induced Pluripotent Stem Cell Generation*. Stem Cells, 2010. 28(8):1315-25.
3. Watanabe, K., M. Ueno, D. Kamiya, A. Nishiyama, M. Matsumura, T. Wataya, J.B. Takahashi, S. Nishikawa, S. Nishikawa, K. Muguruma, and Y. Sasai, *A ROCK inhibitor permits survival of dissociated human embryonic stem cells*. Nat Biotechnol, 2007. 25(6): p. 681-6.
4. Lin, T., R. Ambasudhan, X. Yuan, W. Li, S. Hilcove, R. Abujarour, X. Lin, H.S. Hahm, E. Hao, A. Hayek, and S. Ding, *A chemical platform for improved induction of human iPSCs*. Nat Methods, 2009. 6(11): p. 805-8.

LIMITED USE LICENSE
“For Research Use Only”

This product is to be used for Research Purposes only and for no other purpose. The term “Research Purposes” means scientific research programs directly under the user’s control, which are specifically directed to the purposes of internal research and not for any Commercial Purpose. The purchase and/or use of this product conveys to the user the non-transferable right to use the product and the components of the product in research conducted by the user. The user will not reformulate the Product and shall not sell or otherwise transfer this product, its components, or materials made therefrom to any third party. User shall not use this product or its components for Commercial Purposes until it has obtained a license for this purpose. The term “Commercial Purposes” shall mean any activity by a party for consideration and may include, but is not limited to, use of the product or its components (i) in manufacturing, (ii) to provide a service, information or data to a third party (e.g. fee for sample testing), (iii) for therapeutic, diagnostic or prophylactic purposes, (iv) for resale, whether or not such product or its components are resold for use in research, or (v) in a quality control or quality assurance process for the manufacture of a product for sale. The use of this product by the user constitutes agreement with the terms of this limited use label license for Primorigen Biosciences Inc.”

Primorigen Biosciences Inc. warrants that its products will conform to the standards stated in its product specification sheets in effect at the time of shipment when used as directed. Primorigen Biosciences Inc. will replace, free of charge, any product that does not conform to the specifications, but only when used as directed. This warranty limits the liability of Primorigen Biosciences Inc. to only the replacement of the product. **THIS WARRANTY IS EXCLUSIVE, AND PRIMORIGEN BIOSCIENCES INC. MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The stated express warranties, and the remedy provided for breach thereof, are in lieu of all other liability or obligations of PRIMORIGEN BIOSCIENCES INC. for any damages whatsoever arising out of or in connection with the delivery, use, performance, or the inability to use any of its products. **IN NO EVENT SHALL PRIMORIGEN BIOSCIENCES INC. BE LIABLE UNDER ANY LEGAL THEORY (INCLUDING BUT NOT LIMITED TO CONTRACT, NEGLIGENCE, STRICT LIABILITY IN TORT, OR WARRANTY OF ANY KIND) FOR ANY INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES (INCLUDING BUT NOT LIMITED TO LOST PROFITS).** Without limiting the effect of the preceding sentence, the maximum liability, if any, of Primorigen Biosciences Inc., shall not exceed the purchase price paid by PURCHASER for the product. The use of the product by the user constitutes user's understanding of and agreement with the terms of this limited warranty and limited liability policy.

Primorigen Biosciences • 510 Charmany Drive • Madison, WI 53719

☎ 608.441.8332 • ✉ info@primorigen.com • www.primorigen.com

Trademark and Patent Information: Primorigen®, and Primorigen Biosciences® are registered trademarks of Primorigen Biosciences.

StemAdhere™ is a trademark of Primorigen Biosciences (Patent Pending)